

## BCMA, CD269 anti-human B cell maturation antigen

**Catalog number:** 10100

**名称:** 小鼠抗人B细胞成熟抗原单克隆抗体

**Anti Name:** anti-human B cell maturation antigen

**Source:** Mouse (小鼠)

**Tag:**

**Trait:**

**Concentration:** 1.0mg/ml

ELISA: 1:1,000-1:50,000 (ABTS底物)  
1:5,000-1:100,000 (TMB或OPD底物)  
WB: 1:5000-1:20000

**Immunization:** human B cell细胞成熟抗原重组蛋白

**Research areas:**

**Clone type:** monoclonal (单克隆)

**Kron:** MA807

**Cell:**

**Subtype:** IgG2b( $\kappa$ )

**Weight:**

**Purification:**

**Storage fluid:** 0.15M PBS

**Storage:** -20°C