

## EpCAM/CD326 anti-epithelial cell adhesion molecule/MA8171

Catalog number: 10102

名称: 小鼠抗人上皮细胞黏附分子 (EpCAM/CD326) 单克隆抗体/MA8171

Anti Name: anti-epithelial cell adhesion molecule/MA8171

Source: Mouse (小鼠)

Tag:

Trait:

Concentration: 1.0mg/ml

ELISA: 1:1,000-1:50,000 (ABTS底物)  
1:5,000-1:100,000 (TMB或OPD底物)  
WB: 1:5000-1:20000  
IHC: 1:500-1:2000  
FCM: 1-2μg/反应

Immunization: 人上皮细胞黏附分子重组蛋白

Research areas:

Clone type: monoclonal (单克隆)

Kron: MA8171

Cell:

Subtype: IgG2a(κ)

Weight:

Purification:

Storage fluid: 0.15M PBS

Storage: -20°C