

## EpCAM/CD326 anti-epithelial cell adhesion molecule/MA8172

Catalog number: 10103

**名称:** 小鼠抗人上皮细胞黏附分子 (EpCAM/CD326) 单克隆抗体/MA8172

**Anti Name:** anti-epithelial cell adhesion molecule/MA8172

**Source:** Mouse (小鼠)

**Tag:**

**Trait:**

**Concentration:** 1.0mg/ml

ELISA: 1:1,000-1:50,000 (ABTS底物)  
1:5,000-1:100,000 (TMB或OPD底物)  
WB: 1:5000-1:20000  
IHC: 1:500-1:2000  
FCM: 1-2 $\mu$ g/反应

**Immunization:** 人上皮细胞黏附分子重组蛋白

**Research areas:**

**Clone type:** monoclonal (单克隆)

**Kron:** MA8172

**Cell:**

**Subtype:** IgG1( $\kappa$ )

**Weight:**

**Purification:**

**Storage fluid:** 0.15M PBS

**Storage:** -20°C