

adi anti-human adiponectin

Catalog number: 10105

名称: 小鼠抗人脂联素单克隆抗体

Anti Name: anti-human adiponectin

Source: Mouse(小鼠)

Tag:

Trait:

Concentration: 1.0mg/ml

ELISA: 1:1,000-1:50,000 (ABTS底物)
1:5,000-1:100,000 (TMB或OPD底物)
IP: 1-10 μ g/反应

Immunization: 人脂联素重组蛋白

Research areas:

Clone type: monoclonal(单克隆)

Kron: MA8121

Cell:

Subtype: IgG1(κ)

Weight:

Purification:

Storage fluid: 0.15M PBS

Storage: -20°C