

## RhoE anti-human RhoE

**Catalog number:** 10108

**名称:** 小鼠抗人RhoE单克隆抗体

**Anti Name:** anti-human RhoE

**Source:** Mouse(小鼠)

**Tag:**

**Trait:**

**Concentration:** 1.0mg/ml

ELISA: 1:1,000-1:50,000 (ABTS底物)  
1:5,000-1:100,000 (TMB或OPD底物)  
WB: 1:5000-1:20000  
IHC: 1:500-1:2000  
ICC: 1:500-1:2000

**Immunization:** 人RhoE重组蛋白

**Research areas:**

**Clone type:** monoclonal(单克隆)

**Kron:** MA818

**Cell:**

**Subtype:** IgG1(κ)

**Weight:**

**Purification:**

**Storage fluid:** 0.15M PBS

**Storage:** -20°C