

## ABP anti-human Androgen-binding protein/MA8191

**Catalog number:** 10111

**名称:** 小鼠抗人雄激素结合蛋白ABP单克隆抗体/MA8191

**Anti Name:** anti-human Androgen-binding protein/MA8191

**Source:** Mouse(小鼠)

**Tag:**

**Trait:**

**Concentration:** 1.0mg/ml

ELISA: 1:1,000-1:50,000 (ABTS底物)  
1:5,000-1:100,000 (TMB或OPD底物)  
WB: 1:5000-1:20000  
IHC: 1:500-1:2000

**Immunization:** 人雄激素结合蛋白重组蛋白

**Research areas:**

**Clone type:** monoclonal(单克隆)

**Kron:** MA8191

**Cell:**

**Subtype:** IgG1(κ)

**Weight:**

**Purification:**

**Storage fluid:** 0.15M PBS

**Storage:** -20°C