

FAM9B family with sequence similarity 9 member b

Catalog number: 10114

名称: 小鼠抗人序列相似性家族9成员B单克隆抗体

Anti Name: family with sequence similarity 9 member b

Source: Mouse(小鼠)

Tag:

Trait:

Concentration: 1.0mg/ml

ELISA: 1:1,000-1:50,000 (ABTS底物)
1:5,000-1:100,000 (TMB或OPD底物)
WB: 1:5000-1:20000

Immunization: 抗人序列相似性家族9成员B重组蛋白

Research areas:

Clone type: monoclonal(单克隆)

Kron: MA823

Cell:

Subtype: IgG1(κ)

Weight:

Purification:

Storage fluid: 0.15M PBS

Storage: -20°C